

Phosphate carrier protein, mitochondrial

Catalog No: #AP85233



Package Size: #AP85233-1 50ug #AP85233-2 100ug #AP85233-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphate carrier protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Phosphate transport protein,PTP,Solute carrier family 25 member 3,PHC
Accession No.	Uniprot ID:P12234
Uniprot	P12234
GeneID	282477;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC25A3

Note: This product is for in vitro research use only