

Multidrug resistance-associated protein 9

Catalog No: #AP85234



Package Size: #AP85234-1 50ug #AP85234-2 100ug #AP85234-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Multidrug resistance-associated protein 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATP-binding cassette sub-family C member 12,Mrp9
Accession No.	Uniprot ID:Q80WJ6
Uniprot	Q80WJ6
GeneID	244562;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Abcc12

Note: This product is for in vitro research use only