

Leukocyte-associated immunoglobulin-like receptor 1

Catalog No: #AP85235



Package Size: #AP85235-1 50ug #AP85235-2 100ug #AP85235-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leukocyte-associated immunoglobulin-like receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LAIR-1,CD305
Accession No.	Uniprot ID:Q8BG84
Uniprot	Q8BG84
GeneID	52855;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lair1

Note: This product is for in vitro research use only