

Laminin subunit alpha-1

Catalog No: #AP85236



Package Size: #AP85236-1 50ug #AP85236-2 100ug #AP85236-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Laminin subunit alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Laminin A chain,Laminin-1 subunit alpha,Laminin-3 subunit alpha,S-laminin subunit alpha,S-LAM alpha,LAMA
Accession No.	Uniprot ID:P25391
Uniprot	P25391
GeneID	284217;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LAMA1

Note: This product is for in vitro research use only