

Melanocortin receptor 5

Catalog No: #AP85238



Package Size: #AP85238-1 50ug #AP85238-2 100ug #AP85238-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Melanocortin receptor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MC5-R
Accession No.	Uniprot ID:P56451
Uniprot	P56451
GeneID	507987;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MC5R

Note: This product is for in vitro research use only