

Mitochondrial carnitine/acylcarnitine carrier protein

Catalog No: #AP85239



Package Size: #AP85239-1 50ug #AP85239-2 100ug #AP85239-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitochondrial carnitine/acylcarnitine carrier protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Carnitine/acylcarnitine translocase,CAC,Solute carrier family 25 member 20,Cac,Cact
Accession No.	Uniprot ID:Q9Z2Z6
Uniprot	Q9Z2Z6
GeneID	57279;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc25a20

Note: This product is for in vitro research use only