

Melanin-concentrating hormone receptor 2

Catalog No: #AP85240



Package Size: #AP85240-1 50ug #AP85240-2 100ug #AP85240-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanin-concentrating hormone receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MCH receptor 2,G-protein coupled receptor 145,GPRv17,MCH-2R,MCH2,GPR145,SLT
Accession No.	Uniprot ID:Q969V1
Uniprot	Q969V1
GeneID	84539;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MCHR2

Note: This product is for in vitro research use only