

Ladinin-1

Catalog No: #AP85243



Package Size: #AP85243-1 50ug #AP85243-2 100ug #AP85243-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ladinin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lad-1,Linear IgA disease antigen homolog,LADA
Accession No.	Uniprot ID:P57016
Uniprot	P57016
GeneID	16763;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lad1

Note: This product is for in vitro research use only