

Interstitial collagenase

Catalog No: #AP85244



Package Size: #AP85244-1 50ug #AP85244-2 100ug #AP85244-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interstitial collagenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Fibroblast collagenase,Matrix metalloproteinase-1,MMP-1,CLG
Accession No.	Uniprot ID:P28053
Uniprot	P28053
GeneID	281308;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP1

Note: This product is for in vitro research use only