

Integrin alpha-4

Catalog No: #AP85246



Package Size: #AP85246-1 50ug #AP85246-2 100ug #AP85246-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD49 antigen-like family member D,Integrin alpha-IV,VLA-4 subunit alpha,CD49d,CD49D
Accession No.	Uniprot ID:P13612
Uniprot	P13612
GeneID	3676;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITGA4

Note: This product is for in vitro research use only