

Neural cell adhesion molecule L1

Catalog No: #AP85250



Package Size: #AP85250-1 50ug #AP85250-2 100ug #AP85250-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neural cell adhesion molecule L1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	N-CAM-L1,CD171,Caml1
Accession No.	Uniprot ID:P11627
Uniprot	P11627
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:L1cam

Note: This product is for in vitro research use only