

Lipid phosphate phosphatase-related protein type 4

Catalog No: #AP85253



Package Size: #AP85253-1 50ug #AP85253-2 100ug #AP85253-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lipid phosphate phosphatase-related protein type 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brain-specific phosphatidic acid phosphatase-like protein 1,Plasticity-related gene 1 protein,PRG-1,D3Bwg0562e,Kiaa0455,Prg1
Accession No.	Uniprot ID:Q7TME0
Uniprot	Q7TME0
GeneID	229791;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lppr4

Note: This product is for in vitro research use only