

Lipid phosphate phosphohydrolase 1

Catalog No: #AP85254



Package Size: #AP85254-1 50ug #AP85254-2 100ug #AP85254-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lipid phosphate phosphohydrolase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PAP2-alpha,Phosphatidate phosphohydrolase type 2a,Phosphatidic acid phosphatase 2a,PAP-2a,PAP2a,Lpp1
Accession No.	Uniprot ID:O08564
Uniprot	O08564
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppap2a

Note: This product is for in vitro research use only