

Melanoma-derived growth regulatory protein

Catalog No: #AP85255



Package Size: #AP85255-1 50ug #AP85255-2 100ug #AP85255-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Melanoma-derived growth regulatory protein |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Cartilage-derived retinoic acid-sensitive protein,CD-RAP,Melanoma inhibitory activity protein,Cdrap,Mia1 |
| Accession No. | Uniprot ID:Q61865 |
| Uniprot | Q61865 |
| GeneID | 12587; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mia

Note: This product is for in vitro research use only