

## Mannan-binding lectin serine protease 2

Catalog No: #AP85256



Package Size: #AP85256-1 50ug #AP85256-2 100ug #AP85256-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mannan-binding lectin serine protease 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MBL-associated serine protease 2, Mannose-binding protein-associated serine protease 2, MASP-2
Accession No.	Uniprot ID:O00187
Uniprot	O00187
GeneID	10747;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MASP2

Note: This product is for in vitro research use only