

Muscleblind-like protein 1

Catalog No: #AP85257



Package Size: #AP85257-1 50ug #AP85257-2 100ug #AP85257-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Muscleblind-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Triplet-expansion RNA-binding protein,Exp,Mbnl
Accession No.	Uniprot ID:Q9JKP5
Uniprot	Q9JKP5
GeneID	56758;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mbnl1

Note: This product is for in vitro research use only