

Low-density lipoprotein receptor-related protein 11

Catalog No: #AP85259



Package Size: #AP85259-1 50ug #AP85259-2 100ug #AP85259-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Low-density lipoprotein receptor-related protein 11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LRP-11
Accession No.	Uniprot ID:Q8CB67
Uniprot	Q8CB67
GeneID	237253;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lrp11

Note: This product is for in vitro research use only