

## Lutropin-choriogonadotropic hormone receptor

Catalog No: #AP85262



Package Size: #AP85262-1 50ug #AP85262-2 100ug #AP85262-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lutropin-choriogonadotropic hormone receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LH/CG-R,Luteinizing hormone receptor,LHR,LCGR,LGR2,LHRHR
Accession No.	Uniprot ID:P22888
Uniprot	P22888
GeneID	3973;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LHCGR

Note: This product is for in vitro research use only