

Mitogen-activated protein kinase 8

Catalog No: #AP85265



Package Size: #AP85265-1 50ug #AP85265-2 100ug #AP85265-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP kinase 8, Stress-activated protein kinase JNK1, c-Jun N-terminal kinase 1, Jnk1, Prkm8
Accession No.	Uniprot ID:Q91Y86
Uniprot	Q91Y86
GeneID	26419;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mapk8

Note: This product is for in vitro research use only