

## Matrix metalloproteinase-23

Catalog No: #AP85266



Package Size: #AP85266-1 50ug #AP85266-2 100ug #AP85266-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Matrix metalloproteinase-23
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMP-23,Femalysin,MIFR-1,Matrix metalloproteinase-21,MMP-21,Matrix metalloproteinase-22,MMP-22,MMP21
Accession No.	Uniprot ID:O75900
Uniprot	O75900
GeneID	8510;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MMP23A

Note: This product is for in vitro research use only