

# Molybdenum cofactor sulfurase

Catalog No: #AP85269



Package Size: #AP85269-1 50ug #AP85269-2 100ug #AP85269-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Molybdenum cofactor sulfurase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MCS, Molybdenum cofactor sulfurtransferase, MCSU
Accession No.	Uniprot ID:Q9N0E7
Uniprot	Q9N0E7
GeneID	281226;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MOCOS

Note: This product is for in vitro research use only