

Mitogen-activated protein kinase kinase kinase MLT

Catalog No: #AP85273



Package Size: #AP85273-1 50ug #AP85273-2 100ug #AP85273-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase kinase kinase MLT
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Human cervical cancer suppressor gene 4 protein,HCCS-4,Leucine zipper- and sterile alpha motif-containing kinase,MLK-like mitogen-activated protein triple kinase,Mixed lineage kinase-related kinase,MLK-related kinase,Sterile alpha motif- and leucine zipper-containing kinase AZK,MLTK,HCCS4
Accession No.	Uniprot ID:Q9NYL2
Uniprot	Q9NYL2
GeneID	51776;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ZAK

Note: This product is for in vitro research use only