

Monocarboxylate transporter 1

Catalog No: #AP85277



Package Size: #AP85277-1 50ug #AP85277-2 100ug #AP85277-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Monocarboxylate transporter 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MCT 1,Solute carrier family 16 member 1
Accession No.	Uniprot ID:Q3MHW6
Uniprot	Q3MHW6
GeneID	505775;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC16A1

Note: This product is for in vitro research use only