

Dual specificity mitogen-activated protein kinase kinase 2

Catalog No: #AP85278

Package Size: #AP85278-1 50ug #AP85278-2 100ug #AP85278-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP kinase kinase 2,ERK activator kinase 2,MAPK/ERK kinase 2,MEK 2,Mek2,Mkk2,Prkmk2
Accession No.	Uniprot ID:Q63932
Uniprot	Q63932
GeneID	26396;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Map2k2

Note: This product is for in vitro research use only