

Kallikrein-5

Catalog No: #AP85280



Package Size: #AP85280-1 50ug #AP85280-2 100ug #AP85280-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kallikrein-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kallikrein-like protein 2, KLK-L2, Stratum corneum tryptic enzyme, SCTE, UNQ570/PRO1132
Accession No.	Uniprot ID:Q9Y337
Uniprot	Q9Y337
GeneID	25818;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: KLK5

Note: This product is for in vitro research use only