

Inositol oxygenase

Catalog No: #AP85281



Package Size: #AP85281-1 50ug #AP85281-2 100ug #AP85281-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Inositol oxygenase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Aldehyde reductase-like 6, Myo-inositol oxygenase, MI oxygenase, Renal-specific oxidoreductase, Aldrl6, Rsor |
| Accession No. | Uniprot ID:Q9QXN5 |
| Uniprot | Q9QXN5 |
| GeneID | 56727; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Miox

Note: This product is for in vitro research use only