

Matrix metalloproteinase-16

Catalog No: #AP85282



Package Size: #AP85282-1 50ug #AP85282-2 100ug #AP85282-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MMP-16,Membrane-type matrix metalloproteinase 3,MT-MMP 3,Membrane-type-3 matrix metalloproteinase,MT3-MMP
Accession No.	Uniprot ID:O35548
Uniprot	O35548
GeneID	65205;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmp16

Note: This product is for in vitro research use only