

Lactadherin

Catalog No: #AP85285



Package Size: #AP85285-1 50ug #AP85285-2 100ug #AP85285-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lactadherin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MFGM,Milk fat globule-EGF factor 8,MFG-E8,O-acetyl GD3 ganglioside synthase,AGS,SED1,Ags
Accession No.	Uniprot ID:P70490
Uniprot	P70490
GeneID	25277;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mfge8

Note: This product is for in vitro research use only