

Mitogen-activated protein kinase 8

Catalog No: #AP85286



Package Size: #AP85286-1 50ug #AP85286-2 100ug #AP85286-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP kinase 8,JNK-46,Stress-activated protein kinase 1c,SAPK1c,Stress-activated protein kinase JNK1,c-Jun N-terminal kinase 1,JNK1,PRKM8,SAPK1,SAPK1C
Accession No.	Uniprot ID:P45983
Uniprot	P45983
GeneID	5599;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAPK8

Note: This product is for in vitro research use only