

Microtubule-associated proteins 1A/1B light chain 3A

Catalog No: #AP85287



Package Size: #AP85287-1 50ug #AP85287-2 100ug #AP85287-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Microtubule-associated proteins 1A/1B light chain 3A |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Autophagy-related protein LC3 A, Autophagy-related ubiquitin-like modifier LC3 A, MAP1 light chain 3-like protein 1, MAP1A/MAP1B light chain 3 A, MAP1A/MAP1B LC3 A, Microtubule-associated protein 1 light chain 3 alpha |
| Accession No. | Uniprot ID:Q91VR7 |
| Uniprot | Q91VR7 |
| GeneID | 66734; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Map1lc3a

Note: This product is for in vitro research use only