

Stromelysin-3

Catalog No: #AP85288



Package Size: #AP85288-1 50ug #AP85288-2 100ug #AP85288-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Stromelysin-3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | SL-3,Matrix metalloproteinase-11,MMP-11 |
| Accession No. | Uniprot ID:Q499S5 |
| Uniprot | Q499S5 |
| GeneID | 103694874;25481; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmp11

Note: This product is for in vitro research use only