

Multiple inositol polyphosphate phosphatase 1

Catalog No: #AP85290



Package Size: #AP85290-1 50ug #AP85290-2 100ug #AP85290-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Multiple inositol polyphosphate phosphatase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	2,3-bisphosphoglycerate 3-phosphatase,2,3-BPG phosphatase,Inositol (1,3,4,5)-tetrakisphosphate 3-phosphatase,Ins(1,3,4,5)P(4) 3-phosphatase,Mipp,MNCb-1572
Accession No.	Uniprot ID:Q9Z2L6
Uniprot	Q9Z2L6
GeneID	17330;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Minpp1

Note: This product is for in vitro research use only