

LIM/homeobox protein LMX-1.2

Catalog No: #AP85291



Package Size: #AP85291-1 50ug #AP85291-2 100ug #AP85291-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LIM/homeobox protein LMX-1.2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Homeobox protein LMX-1,LMX,LMX1
Accession No.	Uniprot ID:P53413
Uniprot	P53413
GeneID	396312;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LMX1B

Note: This product is for in vitro research use only