

Muscle, skeletal receptor tyrosine protein kinase

Catalog No: #AP85292



Package Size: #AP85292-1 50ug #AP85292-2 100ug #AP85292-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Muscle, skeletal receptor tyrosine protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Muscle-specific tyrosine protein kinase receptor, MuSK
Accession No.	Uniprot ID:Q8AXY6
Uniprot	Q8AXY6
GeneID	373897;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MUSK

Note: This product is for in vitro research use only