

Myosin light chain 1/3, skeletal muscle isoform

Catalog No: #AP85293



Package Size: #AP85293-1 50ug #AP85293-2 100ug #AP85293-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin light chain 1/3, skeletal muscle isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MLC1/MLC3,Myosin light chain alkali 1/2,Myosin light chain A1/A2
Accession No.	Uniprot ID:P05976
Uniprot	P05976
GeneID	4632;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYL1

Note: This product is for in vitro research use only