

Putative macrophage-stimulating protein MSTP9

Catalog No: #AP85298



Package Size: #AP85298-1 50ug #AP85298-2 100ug #AP85298-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Putative macrophage-stimulating protein MSTP9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Brain rescue factor 1, BRF-1, Hepatocyte growth factor-like protein homolog, D1F15S1A, MSTP9
Accession No.	Uniprot ID: Q2TV78
Uniprot	Q2TV78
GeneID	11223;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MST1P9

Note: This product is for in vitro research use only