

# Inositol oxygenase

Catalog No: #AP85305

Package Size: #AP85305-1 50ug #AP85305-2 100ug #AP85305-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Inositol oxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Aldehyde reductase-like 6, Kidney-specific protein 32, Myo-inositol oxygenase, MI oxygenase, Renal-specific oxidoreductase, ALDRL6, KSP32, RSOR
Accession No.	Uniprot ID: Q9UGB7
Uniprot	Q9UGB7
GeneID	55586;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MIOX

Note: This product is for in vitro research use only