

Proto-oncogene Mas

Catalog No: #AP85306



Package Size: #AP85306-1 50ug #AP85306-2 100ug #AP85306-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Proto-oncogene Mas
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mas,Mas-1
Accession No.	Uniprot ID:P30554
Uniprot	P30554
GeneID	17171;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mas1

Note: This product is for in vitro research use only