

Mitogen-activated protein kinase 1

Catalog No: #AP85312



Package Size: #AP85312-1 50ug #AP85312-2 100ug #AP85312-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MAP kinase 1,ERT1,Extracellular signal-regulated kinase 2,ERK-2,Mitogen-activated protein kinase 2,MAP kinase 2,ERK2,PRKM1
Accession No.	Uniprot ID:P46196
Uniprot	P46196
GeneID	327672;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAPK1

Note: This product is for in vitro research use only