

# Mitogen-activated protein kinase 1

Catalog No: #AP85313



Package Size: #AP85313-1 50ug #AP85313-2 100ug #AP85313-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mitogen-activated protein kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP kinase 1,ERT1,Extracellular signal-regulated kinase 2,ERK-2,MAP kinase isoform p42,p42-MAPK,Mitogen-activated protein kinase 2,MAP kinase 2,Erk2,Mapk,Prkm1
Accession No.	Uniprot ID:P63085
Uniprot	P63085
GeneID	26413;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Mapk1

Note: This product is for in vitro research use only