

Protein lifeguard 1

Catalog No: #AP85316



Package Size: #AP85316-1 50ug #AP85316-2 100ug #AP85316-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein lifeguard 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NMDA receptor glutamate-binding subunit, Putative MAPK-activating protein PM02, NMDARA1
Accession No.	Uniprot ID:Q7Z429
Uniprot	Q7Z429
GeneID	2907;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GRINA

Note: This product is for in vitro research use only