

# LIM homeobox transcription factor 1-beta

Catalog No: #AP85317



Package Size: #AP85317-1 50ug #AP85317-2 100ug #AP85317-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LIM homeobox transcription factor 1-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LIM/homeobox protein 1.2, LMX-1.2, LIM/homeobox protein LMX1B
Accession No.	Uniprot ID:O88609
Uniprot	O88609
GeneID	16917;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Lmx1b

Note: This product is for in vitro research use only