

Lymphatic vessel endothelial hyaluronic acid receptor 1

Catalog No: #AP85318

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85318-1 50ug #AP85318-2 100ug #AP85318-3 1mg

Description

Product Name	Lymphatic vessel endothelial hyaluronic acid receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LYVE-1,Cell surface retention sequence-binding protein 1,CRSBP-1,Extracellular link domain-containing protein 1,Hyaluronic acid receptor,CRSBP1,HAR,XLKD1,UNQ230/PRO263
Accession No.	Uniprot ID:Q9Y5Y7
Uniprot	Q9Y5Y7
GeneID	10894;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LYVE1

Note: This product is for in vitro research use only