

Beta-microseminoprotein

Catalog No: #AP85319



Package Size: #AP85319-1 50ug #AP85319-2 100ug #AP85319-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-microseminoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Prostate secreted seminal plasma protein, Prostate secretory protein of 94 amino acids, PSP-94, TS507, PSP94
Accession No.	Uniprot ID:O02826
Uniprot	O02826
GeneID	396852;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MSMB

Note: This product is for in vitro research use only