

Transcription factor Maf

Catalog No: #AP85323



Package Size: #AP85323-1 50ug #AP85323-2 100ug #AP85323-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor Maf
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Proto-oncogene c-Maf,V-maf musculoaponeurotic fibrosarcoma oncogene homolog,Maf2
Accession No.	Uniprot ID:P54843
Uniprot	P54843
GeneID	17132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Maf

Note: This product is for in vitro research use only