

# Lysophosphatidic acid receptor 1

Catalog No: #AP85326

Package Size: #AP85326-1 50ug #AP85326-2 100ug #AP85326-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lysophosphatidic acid receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	LPA receptor 1,Lysophosphatidic acid receptor Edg-2,Rec1.3,EDG2,LPA1
Accession No.	Uniprot ID:Q28031
Uniprot	Q28031
GeneID	281136;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LPAR1

Note: This product is for in vitro research use only