

Prolow-density lipoprotein receptor-related protein 1

Catalog No: #AP85328



Package Size: #AP85328-1 50ug #AP85328-2 100ug #AP85328-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prolow-density lipoprotein receptor-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LRP-1,Alpha-2-macroglobulin receptor,A2MR,Apolipoprotein E receptor,APOER,CD91,A2MR,APR
Accession No.	Uniprot ID:Q07954
Uniprot	Q07954
GeneID	4035;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LRP1

Note: This product is for in vitro research use only