

## CD74 antibody

Catalog No: #22982

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD74 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 1 and 28 of CD74
Target Name	CD74
Accession No.	Swiss-Prot:P04233Gene ID:972
Uniprot	P04233
GeneID	972;
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


## Application Details

Predicted MW: 18kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:500

## Images




Sample (30 ug of whole cell lysate)

A: Raji

12% SDS PAGE

Primary antibody diluted at 1: 5000



Immunohistochemical analysis of paraffin-embedded SW480 xenograft, using CD74 antibody at 1: 500 dilution.

---

Note: This product is for in vitro research use only