

Protein lifeguard 1

Catalog No: #AP85330



Package Size: #AP85330-1 50ug #AP85330-2 100ug #AP85330-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein lifeguard 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	NMDA receptor glutamate-binding subunit,NMDARA1
Accession No.	Uniprot ID:Q32L53
Uniprot	Q32L53
GeneID	510225;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GRINA

Note: This product is for in vitro research use only