

Keratin, type I cytoskeletal 20

Catalog No: #AP85332



Package Size: #AP85332-1 50ug #AP85332-2 100ug #AP85332-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytokeratin-20,CK-20,Keratin-20,K20
Accession No.	Uniprot ID:Q9D312
Uniprot	Q9D312
GeneID	66809;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Krt20

Note: This product is for in vitro research use only